## **Grids & Specimen Supports**

## **Please Scroll Down**


## **Grid Staining System**

Designed to handle up to 25 grids at one time or as few as a single one. The Matrix has the following features:

- \* Greatly reduces the chance of mechanical damage
- \* Ensures equal staining and rinse times
- \* Has a simple alpha-numeric system for easy identification of each grid.

The plastics used are not acid-resistant and are designed for aqueous solutions only. Note: Do Not use alcohol- based stains; they will damage or destroy the staining vessels. One vessel could, for example, be used for uranyl acetate and the other one for lead citrate.

**G270** Grid staining system **GG270/C** Replacement staing vessels - 1 red and 1 blue



## **Zerostat Anti-Static Gun**

Anti-static generator for neutralising static charges of either polarity on any insulating surface. Use to reduce the risk of grids jumping out of boxes and to eliminate a block charging defect during ultramicrotomy.

G083 Zerostat gun . each



1.21